


# IndiMag 48s

Versatile automated extraction

- Automated, magnetic bead-based extraction for 1 to 48 samples
- Just 30 minutes per run
- Open software for easy protocol customization
- Effective decontamination using UV light
- Compact footprint to fit any lab (14.3 x 16.5 inches/36.4 x 42 cm)
- Optimized plastic use thanks to 3 flexible cartridge/plate formats



## IndiMag Prefilled Cartridges

Simplified nucleic acid extraction

- Save significant hands-on time and plastics with IndiMag prefilled reagent cartridges
- Greater process stability through fewer pipetting steps
- Simply add your samples and lysis buffer and you're ready
- Prefilled to extract 8 to 48 samples per run on IndiMag 48s instruments
- Cartridges also available for other platforms


Get in touch for more information

[www.indical.com](http://www.indical.com)



## Excellent performance you can trust

Figure 1




Maximizing your downstream testing starts with high-quality sample extraction. Our proven IndiMag nucleic acid extraction kits, prefilled reagent cartridges and instruments offer a flexible solution for many starting sample types and pathogens and are suitable for both high- and low-throughput workflows.

**Figure 1.** The IndiMag Pathogen IM48 Cartridge and IndiMag Pathogen Kit deliver comparable excellent performance.

**Figure 2.** IndiMag prefilled cartridges deliver consistent results with flexible batch sizes and walkaway automation.

## Delivering consistent results across platforms\*

Figure 2



\*PRRSV EU serum samples (n = 11)

Request a quote or order online: [shop.indical.com](https://shop.indical.com)

Email: [support@indical.com](mailto:support@indical.com) | [orders@indical.com](mailto:orders@indical.com) | Phone: +49 341 124 54 0 | Fax: +49 341 124 54 60

Email (Americas): [us\\_support@indical.com](mailto:us_support@indical.com) | [us\\_orders@indical.com](mailto:us_orders@indical.com) | Phone (Americas): +1 (407) 766 5128

Connect with us on LinkedIn: [www.linkedin.com/company/indical](https://www.linkedin.com/company/indical)

For up-to-date licensing information and product-specific disclaimers, see the respective handbook or user manual. Regulatory requirements vary by country, products may not be available in your geographic area. Trademarks: INDICAL®, IndiMag®, IndiSpin® (INDICAL BIOSCIENCE GmbH). KingFisher™ (ThermoFisher Scientific). QIAGEN®, QIAcube® (QIAGEN Group). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-IMG48s-IM48KF96-EN-20210930 © 2021 INDICAL BIOSCIENCE

**INDICAL**  
BIOSCIENCE

INDICAL BIOSCIENCE GmbH  
Deutscher Platz 5b  
04103 Leipzig, Germany

INDICAL Inc.  
1317 Edgewater Dr #3722  
Orlando, FL 32804, USA

